

Ready to get started with Dexcom G6?

Welcome to the Dexcom G6 Continuous Glucose Monitoring (CGM) System! The Dexcom G6 app will guide you step-by-step on how to set up your system, insert the sensor and attach the transmitter. It's that simple!

For detailed step-by-step instructions on how to use the Dexcom G6, refer to the [Dexcom G6 Start Here guide](#) in your transmitter box.

COMPONENTS:




Applicator with built-in sensor

Wear for up to 10 days



Transmitter

Lasts up to 90 days



Compatible Smart Device

Check your smart device compatibility online at dexcom.com/compatibility

or



Receiver

If you do not have a compatible smart device*, please contact Dexcom Customer Service to obtain a receiver by calling 1-844-832-1810 (option 1)

GETTING STARTED

1 Download and open the Dexcom G6 app using your compatible smart device*

- For iPhone: Open the App Store, search for Dexcom G6 mmol/L and tap **Get**
- For Android: Open the Google Play Store, search for Dexcom G6 and tap **Install**

2 Follow the on-screen set up instructions



TRAINING RESOURCES

For training videos, handy guides, FAQs and more, scan the QR code or visit

dexcom.com/en-CA/training




NEED MORE ASSISTANCE?

Contact Dexcom CARE to receive personal CGM training and support from Certified Diabetes Educators:

- Setting up your Dexcom G6
- Sensor insertion tips
- Help customizing app settings and features

Call 1-844-832-1810 (Option 4)

Overview of the Dexcom G6 app:



Dexcom Share: Share your glucose data with up to 10 followers

Trend Arrow: Shows what direction your glucose is heading in, and how fast

Current glucose reading

Trend Graph: Shows the last 24 hours of your glucose history data

Settings: Customize your alerts. To change or turn off an alert, go to Settings > Alerts and tap an alert to change it

Smart device sold separately*

Get the most from your Dexcom G6 by using the following apps:



Dexcom Clarity

Use Dexcom Clarity software to see data history and trends. With Clarity you can opt to share data with your healthcare professional and more.



Dexcom Follow[†]

You can let up to 10 people remotely follow your glucose levels. Have your loved ones download the Dexcom Follow App[†], so they can view your glucose data directly from their smart device.

Need more assistance? **Call 1-844-832-1810**

General inquiries: [Select option 1](#)

Insurance questions: [Select option 2](#)

Product replacements & troubleshooting: [Select option 3](#)

New user training & support: [Select option 4](#)

* Smart Device sold separately. For a list of compatible smart devices, please visit dexcom.com/compatibility.

† Separate Follow app and internet connection required. Users should always confirm readings on the Dexcom G6 app or receiver before making treatment decisions.

Dexcom

Dexcom, Dexcom G6, Dexcom Follow, Dexcom Share, and Dexcom Clarity are registered trademarks of Dexcom Inc. in the United States and may be registered in other countries. © 2023 Dexcom Canada, Co. All rights reserved. LBL-1002560 Rev001